



BENCHTOP LOW SPEED CENTRIFUGE BKB-701

LOW SPEED CENTRIFUGE BKB-701

LOW SPEED CENTRIFUGE



Providing multiple self-diagnostic protection functions against door interlock malfunction, over-speed operation, motor overheating, misoperation and communication fault etc. to ensure the user and instrument safety.
With 9 acceleration curves and 10 deceleration curves, it can prevent secondary emulsification or turbidity, and achieve the best centrifugal effect.
Microprocessor control and DC brush-less motor.
RCF value can be calculated automatically. Programmable, save up to 25 user-defined programs.

SPECIFICATIONS

Model	BKB-701
Max. Speed	5500rpm
Max. RCF	5310xg
Max. Capacity	4x500ml
Speed Accuracy	±10r/min
Time Setting Range	1min to 99min59S
Noise	≤65dB(A)
Power Supply	AC220V ±22V 50Hz/60Hz
Power Consumption	750W
Dimension (WxDxH)	470x560x380mm
Packing Dimension (WxDxH)	547x727x422mm
Net Weight	50kg
Alt Name	Low Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Angle Rotor	5500	5310	12x15ml
2	Swing Rotor Yoke	4000	3040	4x500ml
3	Microplate Buckets	4000	3040	4x1x96 well deep 4x2x96 well
4	Swing Rotor Yoke	4000	3040	4x500ml
5	Round buckets	4000	3040	4x500ml
6	Adapter	4000	3040	4x1x250ml
7				4x2x100ml
8				4x4x50ml(conical/round)
9				4x7x20ml
10				4x9x15ml(conical/round)
11				4x11x15ml(conical)
12				4x14x10ml
13				4x19x5/7ml
14				4x36x1.5ml
15	Swing Rotor Yoke	4000	3040	4 rectangular buckets
16	Rectangular buckets			

17	Adapter	4000	3040	4x1x250ml
18				4x3x100ml
19				4x6x50ml(conical/round)
20				4x12x15ml(conical)
21				4x12x15ml(round)
22				4x20x10ml
23				4x31x5/7ml
24				4x56x1.5ml
25	Swing microplate yoke	4000	2200	2x2x96 well
26	Microplate buckets	4000	2200	2x2x96 well
27	Swing Rotor Yoke	5000	4390	4x100ml
28	Brackets			
29	Adapter	5000	4390	4x50ml
30				4x50ml(conical)
31				4x30ml
32				4x20ml
33				4x15ml
34				4x15ml(conical)
35				4x3x5ml blood collection tube
36				4x2x10ml
37				4x4xCumulus tube
38				4x3x1.5ml
39	Swing Rotor Yoke	4200	3100	
40	Brackets	4200	3100	4x2x100ml
41	Adapter	4200	3100	4x2x50ml
42				4x2x50ml(conical)
43				4x2x30ml
44				4x2x20ml
45				4x2x15ml
46				4x2x15ml(conical)
47				4x6x5ml blood collection tube
48				8x2x10ml
49				8x4 cumulus tube
50				8x3x1.5ml
51	Brackets	4200	3100	4x8x15ml
52	Adapter	4200	3100	4x8x10ml
53				4x8x5ml
54	Brackets	4200	3100	4x2x50ml



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com