



BENCHTOP LOW SPEED CENTRIFUGE BKC-1201

LOW SPEED CENTRIFUGE BKC-1201

LOW SPEED CENTRIFUGE LABORATORY



- 1.Adaptability: Versatile fixed-angle rotors catering to various needs.
- 2.Flexible Tube Capacity: Accommodates tubes from 2 ml to 50 ml effortlessly.
- 3.Powerful RCF: Achieves an RCF of over 2,680 x g for precise separations.
- 4.Wide Application Range: From clinical diagnostics to academia and industrial sectors.
- 5.Trusted Performance: Embraced by professionals and researchers globally.

SPECIFICATIONS

Model	BKC-1201
Max. Speed	4000rpm
Max. RCF	2268xg
Max. Capacity	6x50ml
Noise(dBA)	
Dimension(mm)	420x320x270mm
Power supply	AC220V,50HZ,10A, 60W
Speed accuracy	±10rpm
Time range	1min-99minutes
Net weight(without rotor)	18KGS
Alt Name	Low Speed Centrifuge Laboratory

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Max.speed	Max.RCF	Note
1	Angle rotor	24x10/5ml	4000r/min	2268xg	PP/PC tube / Vacuum blood collection tubes
2	Adapter	24x2ml			
3	Adapter	24x1.5ml			
4	Angle rotor	6x50 ml	4000r/min	2268xg	PP/PC tube
5	Adapter	6x15/10ml			
6	Adapter	6x5ml			

APPLICATIONS

Hospitals,Veterinary clinics,Food labs,Universities & colleges, Pharmaceutical labs, Cell culture labs,Infant stations.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com