



-65°C ULTRA LOW TEMPERATURE CHEST FREEZER BKD-402

-60°C ULTRA LOW TEMPERATURE FREEZER BKD-402

-65°C SUPER LOW TEMPERATURE CHEST FREEZER



SPECIFICATIONS

Model	BKD-402
Capacity	50 L
Temperature Range	-25°C ~ -65°C
Climate class	N
Unit dimensions	672*630*785 mm
Interior dimensions	465*350*285 mm
Power supply	220V/50Hz or 115V/60Hz
Rated input	180 W
Power consumption	2.20 kW.h/24h (FRO)
Refrigerants	Mixed CFC free
Compressor	World brand compressor
Thermostat	Brand Electronic controller
Temperature display	Digital display
Evaporator	Wall type copper tube
Condenser	Wall type copper tube
Fan for condenser	German brand, EBM
Insulation thickness	105 mm
Material of interior	SUS 304 stainless steel
Door	Foamed door, top open
Door handle	Recessed handle with lock
Baskets	N/A
Wheels	4 pieces, casters
Net weight	45 kg
Alt Name	-65°C Super Low Temperature Chest Freezer

APPLICATIONS

The -65°C super low chest freezers are used for fresh keeping storage of sashimi such as super frozen tuna, spot prawn, scampi, salmon etc., deep sea foods and high grade foods at super low temperature in the fields of ocean fisheries, Japanese restaurants, seafood markets, hotels & catering industry etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com