



# -65°C ULTRA LOW TEMPERATURE CHEST FREEZER BKD-409

# -60°C ULTRA LOW TEMPERATURE FREEZER BKD-409

-65°C SUPER LOW TEMPERATURE CHEST FREEZER



## SPECIFICATIONS

|                        |                                           |
|------------------------|-------------------------------------------|
| Model                  | BKD-409                                   |
| Capacity               | 460 L                                     |
| Temperature Range      | -25°C ~ -65°C                             |
| Climate class          | N                                         |
| Unit dimensions        | 1400*880*910 mm                           |
| Interior dimensions    | 1225*620*675 mm                           |
| Power supply           | 220 V/50Hz or 115V/60Hz                   |
| Rated input            | 520 W                                     |
| Power consumption      | 6.80 kW.h/24h (FRO)                       |
| Refrigerants           | Mixed CFC free                            |
| Compressor             | World brand compressor                    |
| Temperature controller | Brand Electronic controller               |
| Temperature display    | Digital display                           |
| Evaporator             | Wall type copper tube                     |
| Condenser              | Al finned copper tube                     |
| Fan for condenser      | German brand, EBM                         |
| Insulation thickness   | 88 mm                                     |
| Material of interior   | SUS 304 stainless steel                   |
| Door                   | Foamed door, top open                     |
| Door handle            | Recessed handle with lock                 |
| Baskets                | 1 piece, plastic coated steel wire        |
| Casters                | 4 pieces, casters                         |
| Net weight             | 95 kg                                     |
| Alt Name               | -65°C Super Low Temperature Chest Freezer |

## APPLICATIONS

The -65°C super low chest freezers are used for fresh keeping storage of sashimi such as super frozen tuna, spot prawn, scampi, salmon etc., deep sea foods and high grade foods at super low temperature in the fields of ocean fisheries, Japanese restaurants, seafood markets, hotels & catering industry etc.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)