



-65°C ULTRA LOW TEMPERATURE CHEST FREEZER BKD-413

-60°C ULTRA LOW TEMPERATURE FREEZER BKD-413

-65°C SUPER LOW TEMPERATURE CHEST FREEZER



SPECIFICATIONS

Model	BKD-413
Capacity	1000 L
Temperature Range	-25°C ~ -65°C
Climate class	N
Unit dimensions	2000*1080*910 mm
Interior dimensions	1835*835*690 mm
Power supply	220 V/50Hz or 115V/60Hz
Rated input	850 W
Power consumption	11.50 kW.h/24h (FRO)
Refrigerants	Mixed CFC free
Compressor	World brand compressor
Thermostat	Brand Electronic controller
Temperature display	Digital display
Evaporator	Wall type copper tube
Condenser	Al finned copper tube
Fan for condenser	German brand, EBM
Insulation thickness	85 mm
Material of interior	SUS 304 stainless steel
Door	Foamed door, top open
Door handle	Recessed handle with lock
Baskets	1 piece, plastic coated steel wire
Casters	6 pieces, casters
Net weight	150 kg
Alt Name	-65°C Super Low Temperature Chest Freezer

APPLICATIONS

The -65°C super low chest freezers are used for fresh keeping storage of sashimi such as super frozen tuna, spot prawn, scampi, salmon etc., deep sea foods and high grade foods at super low temperature in the fields of ocean fisheries, Japanese restaurants, seafood markets, hotels & catering industry etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com