



# -65°C ULTRA LOW TEMPERATURE CHEST FREEZER BKD-503

# -60°C ULTRA LOW TEMPERATURE FREEZER BKD-503

-65°C ULTRA LOW TEMPERATURE CHEST FREEZER



## SPECIFICATIONS

|                      |                                           |
|----------------------|-------------------------------------------|
| Model                | BKD-503                                   |
| Capacity             | 300 L                                     |
| Temperature Range    | -25°C ~ -65°C                             |
| Climate class        | N                                         |
| Unit dimensions      | 1690*860*930 mm                           |
| Interior dimensions  | 1010*510*580 mm                           |
| Power supply         | 220V/50Hz or 115V/60Hz                    |
| Rated input          | 400 W                                     |
| Power consumption    | 5.00 kWh/24h (FRO)                        |
| Refrigerants         | Mixed CFC free                            |
| Compressor           | World brand compressor                    |
| Thermostat           | Brand Electronic controller               |
| Temperature display  | Digital display                           |
| Evaporator           | Wall type copper tube                     |
| Condenser            | Wall type copper tube                     |
| Fan for condenser    | German brand, EBM                         |
| Insulation thickness | 125 mm                                    |
| Material of interior | SUS 304 stainless steel                   |
| Door                 | Foamed door, top open                     |
| Door handle          | Recessed handle with lock                 |
| Baskets              | N/A                                       |
| Wheels               | 4 pieces, Foot Master Casters             |
| Net weight           | 150 kg                                    |
| Alt Name             | -65°C Ultra Low Temperature Chest Freezer |

## APPLICATIONS

The -65°C ultra low temperature chest freezers are used for freezing storage of biological samples and products in the fields of universities, institutes, bio-pharmacy, cancer research, vaccine storage and space exploration, many other specialized industries; the freezer also can be used for fresh storage of deep seafood such as tuna, at super low temperature in Japanese restaurants, hotels, high clubs and seafood markets.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)