



BENCHTOP LOW SPEED CENTRIFUGE BKE-2501

LOW SPEED CENTRIFUGE BKE-2501

LOW SPEED TABLE CENTRIFUGE



Intelligent control system
Five groups of programs with humanized design.
Program group shortcut key, can store 99 groups of program , convenient to call at any time.
Automatic electronic door lock.
With instantaneous and continuous centrifugation functions.
Excellent kernel performance
HD large LCD touch display and the parameters can be revised during operation without shutdown.
Microcomputer control ,AC frequency conversion brush-less motor with fast speed and stable operation.
10 gears acceleration /deceleration with linear drive.

SPECIFICATIONS

Model	BKE-2501
Max. Speed	5500rpm
Max. RCF	5071xg
Max. Capacity	4x750ml
Timer	1s~99h59min59s
Revolutions/min	±10r/min
Voltage	AC 220±22V 50Hz 10A or custom
Power	1000W
Noise Level	≤ 62dB(A)
Size	465x580x400 (mm)
Net Weight	62kg
Rotor automatic recognition	Yes
Alt Name	Low Speed Table Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Max Speed(r/min)	Max RCF(xg)
1	Angle rotor(W)	40x15ml(W)	5500rpm	5071xg
2	Angle rotor(W)	6x50ml(W)	5500rpm	3645xg
3	Swing-out rotor	120x2~7ml(Blood collection tube)	4000rpm	3309xg
4		4x750ml(Round bucket)	4000rpm	3309xg
5		4x750ml(Round bucket)/ (with biosafety lid)	4000rpm	3309xg
6		4x750ml Square bucket(used with adapter)	4000rpm	3309xg
7	Swing-out rotor	Auto uncap 104x2-7ml(blood collection tube)	4000rpm	3184xg
8		4x500ml(Round bucket)	4000rpm	3184xg
9		4x500ml(Round bucket) (with biosafetey lid)	4000rpm	3184xg
10		4x500ml Square bucket (used with adapter)	4000rpm	3184xg

11		4x500ml Sharp bucket	4000rpm	3184xg
12	Swing-out rotor	Auto uncap 64x2-7ml(blood collection tube)	4000rpm	3040xg
13		4x250ml(Round bucket)	4000rpm	3040xg
14		4x250ml(Round bucket)/ (with biosafety lid)	4000rpm	3040xg
15		4x250ml Square bucket (used with adapter)	4000rpm	3040xg
16		4x250ml square bucket(with biosafety lid)		
17		4x250ml Sharp bucket(Corning tube)	4000rpm	3040xg
18		6x250ml(Round bucket)	4000rpm	3040xg
19	Swing-out rotor	Multifunctional three-in-one rotor 4x750ml match with below(Bucket/Basket) to use		
20	Hanging bucket	750ml round bucket with biosafety cover include	4000rpm	3425xg
21	Hanging basket	Micro-plate 4x2x96 hole	4000rpm	3425xg
22	Hanging bucket	Auto uncap 120x2-7ml (blood collection tube)	4000rpm	3425xg
23	Hanging bucket	Specialized for 500ml bloodbag (include the bloodbag adapter)	4000rpm	3425xg
24	Swing-out rotor	Micro-plate 2x2x96 hole	4000rpm	2057xg

FEATURES

Intelligent control system

Five groups of programs with humanized design.

Program group shortcut key, can store 99 groups of program , convenient to call at any time.

Automatic electronic door lock.

With instantaneous and continuous centrifugation functions.

Excellent kernel performance

HD large LCD touch display and the parameters can be revised during operation without shutdown.

Microcomputer control ,AC frequency conversion brush-less motor with fast speed and stable operation.

10 gears acceleration /deceleration with linear drive.

Two modes can be set: start timing/arrival revolution timing.

Careful and safe

With rotor automatic identification system .

Over speed,over heat,unbalance protections with auto lid lock.

Steel Structure,centrifuge chamber made of stainless steel with Environmentally friendly anti-corrosion coating treatment.

Intimate detail design

Special noise reduction system.

An intimate tool rack is provided for accompanying tools to guard the instrument at all times.

With an operation SOP guide map, which makes it easier and more standard for laboratory workers.

Equipped with rotor lid that can be quickly locked for easy use.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com