



BENCHTOP LOW SPEED CENTRIFUGE BKE-2601

LOW SPEED CENTRIFUGE BKE-2601

LOW SPEED TABLE CENTRIFUGE



Intelligent control system

Equipped with five groups of user-defined preset programs.

Can store 99 program groups, convenient to call at any time.

Automatic electronic door lock.

One key short-time centrifuge button.

Excellent kernel performance

HD large LCD display and the parameters can be revised during operation without shutdown.

Microcomputer control ,AC frequency conversion brush-less motor with fast speed and stable operation.

10 gear acceleration /deceleration with linear drive.

Two modes can be set: start timing/arrival revolution timing.

SPECIFICATIONS

Model	BKE-2601
Max. Speed	6500rpm
Max. RCF	5091xg
Max. Capacity	4x750ml
Timer	1s~99h59min59s
Revolutions/min	±10r/min
Voltage	AC 220±22V 50Hz 10A or custom
Power	1000W
Noise Level	≤ 62dB(A)
Size	495x625x420 (mm)
Net Weight	70kg
Rotor automatic recognition	Yes
Alt Name	Low Speed Table Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Max Speed(r/min)	Max RCF(xg)
1	Angle rotor(W)	6x50ml(W)	6500rpm	5091xg
2	Angle rotor(W)	40x15ml(W)	5500rpm	5071xg
3	Swing-out rotor	120x2~7ml(Blood collection tube)	4200rpm	3648xg
4		4x750ml(Round bucket)	4200rpm	3648xg
5		4x750ml Square bucket(match below adapter)	4200rpm	3648xg
6	4x750ml square bucket with adapter	4x1x750ml(Low speed flat bottom)	4200rpm	3648xg
7		4x1x500ml(Low speed flat bottom)		
8		4x1x500ml(Conical bottom)		
9		4x1x250ml(Low speed flat bottom)		
10		4x4x100ml(Low speed round bottom tubes)		
11		4x9x50ml(Low speed conical/round bottom universal)		
12		4x19x10ml(Low speed conical/round bottom universal)		

13		4x19x15ml(Low speed conical/round bottom universal)		
14		4x28x10ml Blood collection tubes		
15		4x40x2-7ml Blood collection tubes		
16	Swing-out rotor	Auto uncap 104x2-7ml(blood collection tube)	4200rpm	3510xg
17		4x500ml(Round bucket)	4200rpm	3510xg
18		4x500ml(Round bucket) (with biosafety lid)	4200rpm	3510xg
19		4x500ml Square bucket(used with adapter)	4200rpm	3510xg
20		4x500ml conical bottom(Corning tube)	4200rpm	3510xg
21	Swing-out rotor	Auto uncap 64x2-7ml(blood collection tube)	4200rpm	3352xg
22		4x250ml(Round bucket)	4200rpm	3352xg
23		4x250ml(Round bucket) (with biosafety lid)	4200rpm	3352xg
24		4x250ml Square cup(match with adapter)include adapter: 4x4x50ml conical bottom 4x10x15ml conical bottom	4200rpm	3352xg
25		4x250ml square bucket(with biosafety lid)		
26		4x250ml(Conical bottom)Corning tube	4200rpm	3352xg
27		6x250ml(Round bucket)	4200rpm	3352xg
28	Swing-out rotor body	Multifunctional three-in-one rotor 4x750ml match with below(Bucket/Basket) to use		
29	Hanging bucket	750ml round bucket with biosafety cover	4200rpm	3776xg
30	Hanging basket	Micro-plate 4x2x96 hole	4200rpm	3776xg
31	Hanging bucket	Auto uncap 120x2-7ml (blood collection tube)	4200rpm	3776xg
32	Hanging bucket	Specialized for 500ml bloodbag (include the bloodbag adapter)	4200rpm	3776xg
33	Swing-out rotor	Micro-plate 2x2x96 hole	4200rpm	2267xg

FEATURES

Intelligent control system

Equipped with five groups of user-defined preset programs.
Can store 99 program groups, convenient to call at any time.
Automatic electronic door lock.
One key short-time centrifuge button.

Excellent kernel performance

HD large LCD display and the parameters can be revised during operation without shutdown.
Microcomputer control ,AC frequency conversion brush-less motor with fast speed and stable operation.
10 gear acceleration /deceleration with linear drive.
Two modes can be set: start timing/arrival revolution timing.

Careful and safe

With rotor automatic identification system .
Over speed,over heat,unbalance protections with auto lid lock.
Steel Structure,centrifuge chamber made of stainless steel with Environmentally friendly anti-corrosion coating treatment.

Intimate detail design

An intimate tool rack is provided for accompanying tools to guard the instrument at all times.The front of the fuselage is designed with an operation SOP guide map, which makes it easier and more standard for laboratory personnel to use the instrument.
It is equipped with a rotor cover that can be quickly locked for easy use.
Ergonomic, low opening height, close the centrifuge cover with one press.
Special noise reduction system.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com