



BLOOD CARD CENTRIFUGE BKE-3301

GEL CARD CENTRIFUGE BKE-3301

BLOOD CARD CENTRIFUGE



SPECIFICATIONS

Model	BKE-3301
Max. Speed	4000rpm
Max. RCF	1790xg
Max. Capacity	12
Timer	1min~99min
Revolutions/min	±20r/min
Voltage	AC 220V±22V,50/60Hz,10A
Power	100W
Noise Level	≤ 65dB(A)
Chamber diameter	280mm
Outer Dimensions	385x320x270mm
Package Dimensions	460x390x330 mm
Net weight	15kg
Gross weight	17kg
Alt Name	Blood card centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Speed	RCF
1	Blood card Rotor	12 cards	4000rpm	1790xg



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com