



FLOOR STANDING HIGH SPEED REFRIGERATED CENTRIFUGE BKE-4001

HIGH SPEED REFRIGERATED CENTRIFUGE BKE-4001

FLOOR HIGH SPEED LARGE CAPACITY CENTRIFUGE



This model is small and compact size and able to centrifuge maximum 6x250ml angle rotors. It is the large capacity small sized floor type high speed centrifuge that is not common in the market.

Equipped with both angle rotor for high speed centrifuge and swing-out rotor for low speed centrifuge, multiple functions for one machine.

Microcomputers controlled, AC variable-frequency motor drive, operate stably and quietly

Automatic recognition of rotors and speed limit control, making centrifuging safer. Able to configure centrifugal force, RCF and interchangeable, able to observe anytime

Equipped with short spin press button

SPECIFICATIONS

Model	BKE-4001
Max. Speed	23000rpm
Max. RCF	35500xg
Max. Capacity	4x750ml
timing range	1min~23h59min
Temp Range(freezing)	-20°C~ 40°C
Temp Accuracy(freezing)	±2°C
RPM accuracy	±20r/min
power supply	AC 220±22V 50Hz 15A
total power	1300W
Noise level	≤ 65dB(A)
centrifuge chamber diameter	420mm
dimensions (length x width x height)	510x660x780(mm)
packaging dimensions (length x width x height)	570x720x850(mm)
Net Weight	150kg
Gross Weight	180kg
Alt Name	Floor High Speed Large Capacity Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Speed	RCF
1	Angle Rotor	12x1.5ml	23000rpm	35500xg
2		16x1.5ml	20000rpm	30230xg
3		10/12x5ml	18500rpm	23000xg
4		24x1.5/2.2ml	18000rpm	29650xg
5		12x10ml	17000rpm	29700xg
6		6x50ml/ round bottom	16000rpm	23430xg
7		6x50ml/Sharp bottom	14000rpm	20800xg
8		8x50ml	12000rpm	16000xg
9		4x100ml	13000rpm	18160xg
10		6x100ml	12500rpm	18550xg
11		4x250ml	12000rpm	21630xg

12		4x8x0.2ml PCR Tube	13000rpm	13800xg
13	Swing-out Rotor	96x5/7ml Pipe	4000rpm	3580xg
14		4x500ml	4000rpm	3580xg
15		4x750ml	4000rpm	3580xg
16		6x250ml	4000rpm	3500xg
17	Micro plates rotor	4x2x96 Hole	4000rpm	3550xg

FEATURES

This model is small and compact size and able to centrifuge maximum 6x250ml angle rotors. It is the large capacity small sized floor type high speed centrifuge that is not common in the market.

Equipped with both angle rotor for high speed centrifuge and swing-out rotor for low speed centrifuge, multiple functions for one machine.

Microcomputers controlled, AC variable-frequency motor drive, operate stably and quietly

Automatic recognition of rotors and speed limit control, making centrifuging safer.

Able to configure centrifugal force, RCF and interchangeable, able to observe anytime

Equipped with short spin press button

Second display for countdown less than one minute

Adopted electronic door lock, inner chamber protected by steel materials

Suitable for radioactive immunology, biochemistry, pharmaceuticals, isolation and purification of blood samples.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com