



BENCHTOP HIGH SPEED CENTRIFUGE BKE-4301

HIGH SPEED CENTRIFUGE BKE-4301

BENCHTOP HIGH SPEED CENTRIFUGE



High and low speed centrifuge, equipped with multiple types of rotors and stands, superior compatibility and multifunction

Controlled by computers, AC frequency variable motor drive, able to operate stably and quietly

Multi-color LCD display, user-friendly, clearer and more direct display.

Able to configure centrifugal force, RCF and interchangeable, able to observe anytime

10 types of acceleration and deceleration control, the 9th control able to achieve free stopping time longer than 540s, able to satisfy the requirement of some special samples.

Second display for countdown less than one minute

Adopted electronic door lock, inner chamber protected by steel materials

SPECIFICATIONS

Model	BKE-4301
Max. Speed	21000rpm
Max. RCF	30642xg
Max. Capacity	4x750ml
Timer	1min~99min
Revolutions/min	±20r/min
Voltage	AC 220±22V 50Hz 10A
Power	750W
Chamber diameter	≤ 65dB(A)
Noise Level	420mm
Outer Dimensions	560x460x400mm
Outer packing Dimensions	640x530x470mm
Net weight	45kg
Gross weight	50kg
Alt Name	Benchtop High Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Speed	RCF
1	Angle Rotor	12x1.5/2.2ml	21000rpm	30642xg
2		16x1.5/2.2ml	20000rpm	30230xg
3		18x0.5ml	18000rpm	20200xg
4		10/12x5ml	18000rpm	20800xg
5		24x1.5/2.2ml	16000rpm	23200xg
6		12x10ml	15000rpm	23120xg
7		12x15ml /Sharp bottom	10000rpm	9690xg

8		6x50ml/Round bottom	14000rpm	19885xg
9		6x50ml/Sharp bottom	13000rpm	18065xg
10		2x8x0.2ml PCR Tube	13500rpm	13880xg
11		4x8x0.2ml PCR Tube	13000rpm	13880xg
12		8x50ml	12000rpm	16000xg
13		4x100ml	11000rpm	13000xg
14		6x100ml	10000rpm	11960xg
15	Swing-out Rotor	4x750ml	4200rpm	3950xg
16		4x500ml	4200rpm	3950xg
17		96x5/7ml	4200rpm	3950xg
18		4x250ml	4200rpm	3850xg
19		6x250ml	4200rpm	3850xg
20		4x50/100ml	5000rpm	4390xg
21		8x50/100ml	4200rpm	3100xg
22		16x10/15ml	4200rpm	3100xg
23		24x10/15ml	4200rpm	3100xg
24		32x10/15ml	4200rpm	3100xg
25	Elisa plate rotor	4x2x96 Hole	4200rpm	3910xg
26		2x2x96 Hole	4200rpm	2530xg
27		2x2x48 Hole	4200rpm	1900xg

FEATURES

High and low speed centrifuge, equipped with multiple types of rotors and stands, superior compatibility and multifunction

Controlled by computers, AC frequency variable motor drive, able to operate stably and quietly

Multi-color LCD display, user-friendly, clearer and more direct display. .

Able to configure centrifugal force, RCF and interchangeable, able to observe anytime

10 types of acceleration and deceleration control, the 9th control able to achieve free stopping time longer than 540s, able to satisfy the requirement of some special samples.

Second display for countdown less than one minute

Adopted electronic door lock, inner chamber protected by steel materials

Suitable for radioactive immunology, biochemistry, pharmaceuticals, isolation and purification of blood samples.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com