



CELL WASHING CENTRIFUGE BKE-4501

CELL WASHING CENTRIFUGE BKE-4501

TABLE CELL WASHING CENTRIFUGE



Small size, great space saver for the lab.
Steel structure, centrifuge chamber made of stainless steel.
AC frequency variable motor drives, stably and quietly during operation.
Parameters can be revised during operation without stopping.
Automatic calculates and displays RCF at the same time.

SPECIFICATIONS

Model	BKE-4501
Max. Speed	3000rpm
Max. RCF	890xg
Max. Capacity	12 cards
Timer	1 min ~ 99 min
Revolutions/min	±10 r/min
Voltage	AC 220 ± 22 V, 50/60 Hz, 10 A
Power	100 W
Noise Level	≤ 60 dB(A)
Chamber Diameter	280 mm
Dimensions (L x W x H)	385 x 320 x 270 mm
Packaging Dimensions	460 x 390 x 330 mm
Net Weight	15 kg
Alt Name	Table Cell Washing Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Speed	RCF
1	Cell Washing Rotor	12x(0.25-1)ml	4700rpm	2000xg
2		12x7ml	2976rpm	1000xg



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com