



BENCHTOP LOW SPEED CENTRIFUGE BKE-4801

LOW SPEED CENTRIFUGE BKE-4801

LOW SPEED CENTRIFUGE



Small footprint; great space saver for the lab
Steel structure, centrifuge chamber made of stainless steel.
Maintenance free brushless DC motor drive to enable stable and quiet operation
Multi-color LCD display, user-friendly, clearer and more direct display. .
Able to change parameters anytime during operation without stopping the machine
Automatic calculation and simultaneous display of RCF/RPM.
10 types of accelerating/decelerating control, 40 groups of programming space, user able to program and call program freely
With electronic door lock, enhanced safety
This model is widely used in laboratories in hospitals and institutions due to its super high performance cost ratio and reliable quality

SPECIFICATIONS

Model	BKE-4801
Max. Speed	4000rpm
Max. RCF	2250xg
Max. Capacity	4/6x50ml
Timer	1min~99min
Revolutions/min	±20r/min
Voltage	AC 220V±22V,50/60Hz,10A
Power	100W
Noise Level	≤ 65dB(A)
Chamber diameter	280mm
Outer Dimensions	385x320x270mm
Outer packing Dimensions	460x390x330 mm
Net weight	15kg
Gross weight	17kg
Alt Name	Low speed centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Max. Speed	Max. RCF
1	Angle Rotor	4/6x50ml	4000rpm	2250xg
2		12x10 ml	4000rpm	2250xg
3		12x15 ml	4000rpm	2250xg
4		18x10ml	4000rpm	2250xg
5		24x10ml	4000rpm	2250xg
6		12x20ml	4000rpm	2250xg



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com