



BENCHTOP LOW SPEED CENTRIFUGE BKE-5001

LOW SPEED CENTRIFUGE BKE-5001

LOW SPEED CENTRIFUGE



Compact size; great space saver for the lab

Steel structure, centrifuge chamber made of stainless steel.

Controlled by microcomputers, AC frequency variable motor drive, able to operate stable and quietly.

Multi-color LCD display, user-friendly, clearer and more direct display. .

Able to change parameters anytime during operation without stopping the machine.

Automatic calculation and simultaneous display of RCF/RPM.

10 types of accelerating/decelerating control, 40 groups of programming space, user able to program and call program freely.

With electronic door lock, enhanced safety.

This model is widely used in laboratories in hospitals and institutions due to its super high performance cost ratio and reliable quality.

SPECIFICATIONS

Model	BKE-5001
Max. Speed	5000rpm
Max. RCF	4020xg
Max. Capacity	4x100ml
Timer	1min~99min
Revolutions/min	±20r/min
Voltage	AC 220V±22V,50/60Hz,10A
Power	150W
Noise Level	≤ 65dB(A)
Chamber diameter	320mm
Outer Dimensions	460x360x325mm
Outer packing Dimensions	545x430x395mm
Net weight	35kg
Gross weight	40kg
Alt Name	Low speed centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Max. Speed	Max. RCF
1	Angle Rotor(W)	4/6x50ml	5000rpm	4020xg
2		12x10ml	5000rpm	4020xg
3		12x15ml	5000rpm	4020xg
4	Angle Rotor	12x10 ml	4000rpm	2250xg
5		12x15 ml	4000rpm	2250xg
6		12x20 ml	4000rpm	2250xg
7		18x10ml	4000rpm	2250xg
8		24x10ml	4000rpm	2250xg
9	Swing-out Rotor(Y)	12/24x5ml	4000rpm	2250xg
10		12/24x10ml	4000rpm	2810xg
11		4/6x50ml	4000rpm	2810xg
12		8/16x10/15ml	4000rpm	2810xg
13		4x100ml	4000rpm	2810xg
14	Swing-out Rotor	16x15ml+4x50ml	4000rpm	2810xg

15		4x50ml	4000rpm	2810xg
----	--	--------	---------	--------



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
 Email: info@biolabscientific.com | Website: www.biolabscientific.com