



STATIC PASS BOX BKF-1403

PASS BOX BKF-1403

GMP STANDARD PASS BOX



Fully made of stainless steel 304
Front windows are made of tempered glass
Equipped with UV lamp in the work area
With micro-computer interlock to ensure two doors won't open at the same time
Dynamic pass box equipped with the fan and HEPA filter

SPECIFICATIONS

Model	BKF-1403
Internal size mm (WxDxH)	600x600x600
External size mm (WxDxH)	785x666x720
Alt Name	GMP Standard Pass Box



APPLICATIONS

The pass box is an auxiliary equipment for clean rooms, mainly used for the transfer of small items between clean and non clean areas, in order to reduce the number of open doors in the clean room and minimize pollution to the clean room.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com