



STATIC PASS BOX BKF-1407

PASS BOX BKF-1407

GMP STANDARD PASS BOX



Fully made of stainless steel 304
Front windows are made of tempered glass
Equipped with UV lamp in the work area
With micro-computer interlock to ensure two doors won't open at the same time
Dynamic pass box equipped with the fan and HEPA filter

SPECIFICATIONS

Model	BKF-1407
Internal size mm (WxDxH)	1000x1000x1000
External size mm (WxDxH)	1185x1066x1120
Alt Name	GMP Standard Pass Box



APPLICATIONS

The pass box is an auxiliary equipment for clean rooms, mainly used for the transfer of small items between clean and non clean areas, in order to reduce the number of open doors in the clean room and minimize pollution to the clean room.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com