



## DYNAMIC PASS BOX BKF-1413

## PASS BOX BKF-1413

GMP STANDARD PASS BOX



Fully made of stainless steel 304  
Front windows are made of tempered glass  
Equipped with UV lamp in the work area  
With micro-computer interlock to ensure two doors won't open at the same time  
Dynamic pass box equipped with the fan and HEPA filter

## SPECIFICATIONS

|                          |                       |
|--------------------------|-----------------------|
| Model                    | BKF-1413              |
| Internal size mm (WxDxH) | 600x600x600           |
| External size mm (WxDxH) | 785x666x1100          |
| Alt Name                 | GMP Standard Pass Box |



## APPLICATIONS

The pass box is an auxiliary equipment for clean rooms, mainly used for the transfer of small items between clean and non clean areas, in order to reduce the number of open doors in the clean room and minimize pollution to the clean room.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)