



AIR SAMPLER BKF-2401

AIR SAMPLER BKF-2401

BM-100 PLANKTONIC BACTERIA SAMPLER CALIBRATOR



SPECIFICATIONS

Model	BKF-2401
High precision	±2.5%
Flow range	0-200L/min
Appearance size	diameter 110mm, height 87mm
Operation Method	One-touch operation,you can use it when you turn it on
lithium battery	7.4V 3000AH
charger	8.4V 1A
Response quickly	data 1S is updated once
Alt Name	Planktonic bacteria sampler calibrator

APPLICATIONS

Air bacteria sampler calibratorMainly used to calibrate all microporous plankton samplers.The sampler flow can be measured quickly and accurately without increasing the load on the sampler, causing large deviations in the flow before and after the test. The accuracy of the data can be comparable to that of the DA-100NT.

Executive standard for : JJG198-1994.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com