



AIR SAMPLER BKF-2701

AIR SAMPLER BKF-2701

FKC-V ON-LINE AIR SAMPLER



SPECIFICATIONS

Model	BKF-2701
Volume	100 L/min
Sampling volume	Settable
Petri dish specifications	(ϕ 90 x 15) mm
Power supply	DC24V, 5A
Sampling head	Micropores, alumina
Overall dimensions	224130245 (mm) LWH
Net weight	2.7g
Power dissipation	<120 W
Alt Name	FKC-V on-line Air sampler

APPLICATIONS

The sampler consists of two parts: the fan box and the sampling head. The porous sampling head is made of stainless steel 316L or alumina, and there are 397 0.7mm microholes on it.

Optional bottom bracket to meet different sampling heights.

High performance fans meet higher impact speeds, ensuring the capture efficiency of smaller particle size microorganisms.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com