



AEROSOL GENERATOR BKF-3901

AEROSOL GENERATOR BKF-3901

ZJSJ-012 AEROSOL GENERATOR



SPECIFICATIONS

Model	BKF-3901
Aerosol output volume flow rate	1-2.5m ³ /h
Aerosol production rate	0.5-2x10 ¹³ P/h
Atomization pressure	≥10Kpa - 80Kpa adjustable
Particle size distribution of aerosol	0.02-2um
Aerosol type	DEHS cold polydispersity, PAO/DOP cold polydispersity
Size HxWxL	H200xW250xL370 mm
Weight(kg)	3kg
External air source (compressed air)	Clean and stable gas supply less than 0.7Mpa
Alt Name	aerosol generator

APPLICATIONS

aerosol generator is a portable cold generated aerosol generator developed and produced by our company. This aerosol generator is suitable for air flow rates of ≤50000m³/h and can produce a maximum output of 0.5-2 x 10¹³ P/h polydispersible aerosol. It can meet the requirements of HEPA filter testing for all small clean units to large clean room air supply systems. Produces a continuous and stable aerosol particle, which is convenient for adding solution.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com