



LASER DUST PARTICLE COUNTER BKF-4201

LASER DUST PARTICLE COUNTER BKF-4201

ZJSJ-010 PARTICLE DILUTER



SPECIFICATIONS

Model	BKF-4201
Dilution ratio	60-700 times adjustable @28.3L/min; 10-80 times adjustable @2.83L/min, other magnification can be customized.
Counting flowmeter	28.3L/min
Dilution ratio control mode	Flow regulation mode
Alt Name	Particle diluter

APPLICATIONS

The particle diluter is a portable diluter manufactured by the company. The dilution ratio of this diluter can be adjusted to meet the needs of different customers in different environments, covering a range of HEPA filter testing from small clean units to large clean room air supply systems. The diluter can form a stable dilution ratio and effectively protect the counter.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com