



VISIBLE SPECTROPHOTOMETER BKG-201

VISIBLE SPECTROPHOTOMETER BKG-201

VISIBLE SPECTROPHOTOMETER



SPECIFICATIONS

Model	BKG-201
Wavelength Range	320-1000nm
Spectral Bandwidth	6nm
Wavelength Accuracy	±2nm
Light Source Lamp	Halogen-Tungsten Lamp, 20 W / 12 V
Wavelength Reproducibility	1 nm
Photometric Accuracy	±2.0% (T) (SRM930D)
Photometric Reproducibility	0.5% (T)
Stray Light	≤0.5% (T)
Noise	±0.5% (T)
Value Range - TRANS	0-199.9%
Value Range - ABS	-0.3-2.999
Value Range - FACT	1-9999
Value Range - CONC	0-9999
Power Requirement	100-240 V, 50/60 Hz
Interface	RS232 Serial
Alt Name	Visible Spectrophotometer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com