





STEEL WOOD H FRAME BENCH BKJ-303


STEEL WOOD BENCH BKJ-303

STEEL-WOOD LAB SIDE BENCH



SPECIFICATIONS

Model	BKJ-303
Size (W*D*H)	1000*750*850 mm
Structure	H-frame steel square tube with all-wood base cabinet structure
C-Frame	40*60 mm x 1.2 mm steel square tube
Cabinet & Drawer	15 mm medium-density fiberboard with 1 mm PVC edge sealing
Handle	Embedded aluminum alloy handle with epoxy coating
Table Top	12.7 mm thick black phenolic resin top
Drawer Slides	Three-section silent drawer slides
Hinges	Soft-close hinges
Configuration	Single cabinet with workspace; cabinet can be on the left or right side
Alt Name	Steel-wood lab side bench



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com