



PP BENCH BKJ-401

PP BENCH BKJ-401

ALL PP LAB SIDE BENCH



SPECIFICATIONS

Model	BKJ-401
Size (W*D*H)	1500*750*850mm
Structure	All PP, floor-standing structure
Cabinet and Drawer	8 mm PP (polypropylene) sheet, welded using rods of the same color and material
Table Top	12.7 mm thick black phenolic resin top
Drawer Slides	PP drawer slides
Hinges	PP hinges
Handle	C-type or all-PP straight pull handle
Water Tap	_____
Water Sink	_____
Configuration	Two cabinets - left and right cabinet
Alt Name	All PP lab side bench



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com