



CLASS B AUTOCLAVE BKK-1102

CLASS B AUTOCLAVE BKK-1102

PORTABLE PRESSURE STEAM STERILIZER



Standard Class B with three times vacuum and drying, the remaining temperature of the sterilized instrument is less than 0.2%.

LCD display working status, touch keys.

The vacuum data can reach -0.8 bar, it fits various instruments, including package, unpackaged, solid, hollow, multi-apertures, and inset pipes. It makes sure the instrument can be sterilized completely.

Fully computer-controlled LCD display bright show. The interface is a modularized panel for easy operation.

It is installed with BOWIE&DICK which measures the penetration of water steam.

It has a vacuum test which can ensure the ability of vacuum testing.

Fast and independent steam generator.

Built-in mini printer which can record the sterilizing information.

SPECIFICATIONS

Model	BKK-1102
Capacity	24L
Sterilization Temperature	134°C
Working Pressure	0.22 Mpa
Features	
Sterilizing chamber volume	24L Φ250x470mm
Temperature adjust range	105-134°C
Timer range	0-99min
Heat average	±1°C
Power	AC220V 50Hz/2KW
Dimension(mm)	480x680x455
Exterior package dimension	580x790x580
G.W/N.W	57kg/53kg
Alt Name	Portable Pressure Steam Sterilizer

FEATURES

Standard Class B with three times vacuum and drying, the remaining temperature of the sterilized instrument is less than 0.2%.

LCD display working status, touch keys.

The vacuum data can reach -0.8 bar, it fits various instruments, including package, unpackaged, solid, hollow, multi-apertures, and inset pipes. It makes sure the instrument can be sterilized completely.

Fully computer-controlled LCD display bright show. The interface is a modularized panel for easy operation.

It is installed with BOWIE&DICK which measures the penetration of water steam.

It has a vacuum test which can ensure the ability of vacuum testing.

Fast and independent steam generator.

Built-in mini printer which can record the sterilizing information.

Equipped with Safety valve, Safety lock door system, Pressure lock system, Pressure or temperature overload protection, and Alarm system.

It inactivates the HIV, HBV, BSE, and bacillus efficiently.

This product is designed for open water tanks, water purification, and sewage tank separation, convenient and intuitive.

APPLICATIONS

The table type sterilizers are safely auto-controlled, which are designed for the dental clinics, hospitals, school lab, scientific research institutions and laboratory. And mostly suitable for sterilizing surgical, dental and eye instruments, glassware medicine, culture medium and biological dressing, food, etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com