



BENCHTOP AUTOCLAVE BKK-401

BENCHTOP AUTOCLAVE BKK-401

BENCHTOP STEAM STERILIZER



For 4~6 minutes rapidly sterilize.
Digital display of working status, touch type key.
With 3 fixed cycles of the sterilizing program and adjustable for the user.
The whole process of water adding, temperature rising, sterilizing, drying steam discharging is controlled automatically.
Steam-water inner circulation system: no steam discharge and the environment for sterilizing will be clean and dry.
Safe protection of water lacking.
Exhaust the cool air automatically.
Door safety lock system.
With four stainless steel sterilizing plates.
The chamber of the sterilizer is made of stainless steel.

SPECIFICATIONS

Model	BKK-401
Capacity	35L
Sterilization Temperature	134°C
Working Pressure	0.22 Mpa
Features	With drying function
Sterilizing chamber volume	35L Φ 300x500mm
Temperature adjust range	105-134°C
Timer range	0-99min
Heat average	$\leq \pm 1^\circ\text{C}$
Power	2.7KW/AC220V 50Hz
Sterilizing plate	400x200x30(4piece)
Dimension	730x550x510
Exterior package dimension	780x670x650
G.W/N.W	92kg/78Kg
Alt Name	Benchtop Steam Sterilizer

FEATURES

For 4~6 minutes rapidly sterilize.
Digital display of working status, touch type key.
With 3 fixed cycles of the sterilizing program and adjustable for the user.
The whole process of water adding, temperature rising, sterilizing, drying steam discharging is controlled automatically.
Steam-water inner circulation system: no steam discharge and the environment for sterilizing will be clean and dry.
Safe protection of water lacking.
Exhaust the cool air automatically.
Door safety lock system.
With four stainless steel sterilizing plates.
The chamber of the sterilizer is made of stainless steel.
Automatically shut off with beep reminding after sterilization.
With a drying function.

APPLICATIONS

The table type sterilizers are safely auto-controlled, which are designed for the dental clinics, hospitals, school lab, scientific research institutions and laboratory. And mostly suitable for sterilizing surgical, dental and eye instruments, glassware medicine, culture medium and biological dressing, food, etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com