



CLASS B AUTOCLAVE BKK-702

CLASS B AUTOCLAVE BKK-702

PULSE VACUUM STEAM STERILIZER



Standard Class B with three times vacuum and drying, the remaining temperature of the sterilized instrument is less than 0.2%.

LCD display working status, touch keys.

The vacuum data can reach -0.8 bar, it fits various instruments, including package, unpackaged, solid, hollow, multi-apertures, and inset pipes. It makes sure the instrument can be sterilized completely.

Fully computer-controlled LCD display bright show. The interface is a modularized panel for easy operation.

It is installed with BOWIE & DICK which measures the penetration of water steam.

It has a vacuum test which can ensure the ability of vacuum testing.

Fast and independent steam generator.

Built-in mini printer which can record the sterilizing information.

SPECIFICATIONS

Model	BKK-702
Capacity	50L
Sterilization Temperature	134°C
Working Pressure	0.22 Mpa
Features	
Sterilizing chamber volume	50L Φ 340x550mm
Temperature adjust range	105-134°C
Timer range	0-99min
Heat average	$\pm 1^\circ\text{C}$
Power	AC220V 50Hz/2.7KW
Dimension(mm)	600x850x560
Exterior package dimension	690x900x700
G.W/N.W	125kg/110kg
Alt Name	Pulse Vacuum Steam Sterilizer

FEATURES

Standard Class B with three times vacuum and drying, the remaining temperature of the sterilized instrument is less than 0.2%.

LCD display working status, touch keys.

The vacuum data can reach -0.8 bar, it fits various instruments, including package, unpackaged, solid, hollow, multi-apertures, and inset pipes. It makes sure the instrument can be sterilized completely.

Fully computer-controlled LCD display bright show. The interface is a modularized panel for easy operation.

It is installed with BOWIE & DICK which measures the penetration of water steam.

It has a vacuum test which can ensure the ability of vacuum testing.

Fast and independent steam generator.

Built-in mini printer which can record the sterilizing information.

Equipped with Safety valve, Safety lock door system, Pressure lock system, Pressure or temperature overload protection, and Alarm system.

It inactivates the HIV, HBV, BSE, and bacillus efficiently.

APPLICATIONS

The table type sterilizers are safely auto-controlled, which are designed for the dental clinics, hospitals, school lab, scientific research institutions and laboratory. And mostly suitable for sterilizing surgical, dental and eye instruments, glassware medicine, culture medium and biological dressing, food, etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com