



## LABORATORY VERTICAL AUTOCLAVE BKK-802

# LABORATORY VERTICAL AUTOCLAVE BKK-802

## VERTICAL PULSE VACUUM STEAM STERILIZER



7" Capacitive touch screen, fully controlled automatically, easy to operate. Raising the door, which is more convenient to place and retrieve the items from the autoclave.

The whole procedure from water inlet, pulse vacuum, heating, sterilization, drainage, exhausting to vacuum drying, which is all automatic by one key operation.

Pulse vacuum, the limit is -0.080MPa, which ensures the sterilization more completely.

Separated steam generator, which is more energy-saving compared with normal heating and prevents the sterilizing items from being oxidized.

There are multiple programs loaded in the autoclaves for different requirements of different items.

With the PT-TT test port, it is easy to run the testing.

## SPECIFICATIONS

|                              |                                                       |
|------------------------------|-------------------------------------------------------|
| Model                        | BKK-802                                               |
| Capacity                     | 80L                                                   |
| Sterilization Temperature    | 134°C                                                 |
| Working Pressure             | 0.22 MPa                                              |
| Features                     |                                                       |
| Chamber Volume               | 80LΦ396 x 650mm                                       |
| Height after climb (mm)      | 1590                                                  |
| Heat average                 | ± 1 °C                                                |
| Timer range                  | 0 ~ 99min                                             |
| Temperature adjustment range | 105 ~ 134 °C                                          |
| Drive                        | Steam generator: 3kW Dry blanket: 2.4 kW/AC220V .50Hz |
| Overall dimension (mm)       | 686 x 556 x 1035                                      |
| Transport dimension (mm)     | 786 x 656 x 1135                                      |
| G.win.w                      | 230 kg/200 kg                                         |
| Alt Name                     | Vertical Pulse Vacuum Steam Sterilizer                |

## FEATURES

7" Capacitive touch screen, fully controlled automatically, easy to operate.

Raising the door, which is more convenient to place and retrieve the items from the autoclave.

The whole procedure from water inlet, pulse vacuum, heating, sterilization, drainage, exhausting to vacuum drying, which is all automatic by one key operation.

Pulse vacuum, the limit is -0.080MPa, which ensures the sterilization more completely.

Separated steam generator, which is more energy-saving compared with normal heating and prevents the sterilizing items from being oxidized.

There are multiple programs loaded in the autoclaves for different requirements of different items.

With the PT-TT test port, it is easy to run the testing.

Equipped with Germany Sartorius Sterile Ventilation Filter Midisart 2000 0.20µm, which ensures the air is sterile and prevents pollution.

Steam-Water Inner Circulation System: No steam emission and the environment for sterilization will be clean and dry.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)