



LABORATORY CHAIR BKL-1003


LABORATORY CHAIR BKL-1003

LABORATORY CHAIR



SPECIFICATIONS

Model	BKL-1003
Material	100% Polyurethane
Seat Size(W*D)	415*390mm
Back Size(W*H)	380*250mm
Function	Backrest forward and backward adjustment/Seat up and down adjustment
Footrest	Φ450mm/ High quality, 201 stainless steel
Seat Height	525-715mm
Base	Φ640mm / High-quality cold-rolled steel ,chrome finished
Wheels	PA Nylon, Fixed and Moveable optional
Color	Black
Alt Name	Laboratory Chair



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com