



LABORATORY CHAIR BKL-1203

LABORATORY CHAIR BKL-1203

LABORATORY CHAIR



SPECIFICATIONS

Model	BKL-1203
Material	100% Polyurethane
Seat Size(W*D)	355*310mm
Back Size(W*H)	135mm
Function	Seat up and down adjustment
Footrest	Φ350mm/ High quality, 201stainless steel
Seat Height	530-780mm
Base	Φ480mm / High-quality cold-rolled steel ,chrome finished
Wheels	PA Nylon, Fixed and Moveable optional
Color	Black
Alt Name	Laboratory Chair



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com