



LABORATORY CHAIR BKL-1704


LABORATORY CHAIR BKL-1704

LABORATORY CHAIR



SPECIFICATIONS

Model	BKL-1704
Material	100% Polyurethane
Seat Size(W*D)	470*405mm
Function	Seat up and down adjustment
Footrest	Φ400mm/ High quality, 201stainless steel
Seat Height	595-845mm
Base	Φ640mm / Die casting aluminum alloy ,polishing finished
Wheels	PA Nylon, Fixed and Moveable optional
Color	Black
Alt Name	Laboratory Chair



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com