



## LABORATORY CHAIR BKL-1804


# LABORATORY CHAIR BKL-1804

## LABORATORY CHAIR



## SPECIFICATIONS

|                |                                                         |
|----------------|---------------------------------------------------------|
| Model          | BKL-1804                                                |
| Material       | 100% Polyurethane                                       |
| Seat Size(W*D) | 420*340mm                                               |
| Function       | Seat up and down adjustment                             |
| Footrest       | Φ400mm/ High quality ,201stainless steel                |
| Seat Height    | 595-845mm                                               |
| Base           | Φ640mm / Die casting aluminum alloy ,polishing finished |
| Wheels         | PA Nylon, Fixed and Moveable optional                   |
| Color          | Black                                                   |
| Alt Name       | Laboratory Chair                                        |



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)