



DUCTLESS FUME HOOD BKL-201

DUCTLESS FUME HOOD BKL-201

BENCH-TOP FUME HOOD



SPECIFICATIONS

Model	BKL-201
Internal Size (W*D*H)	960*680*1150mm
External Size (W*D*H)	1200*850*1500mm
Exhaust air volume	864-1296m³/h
Air Velocity	0.4-0.6 m/s
Material	1.0mm full-thickness high-quality cold-roll steel sheet
Lighting	LED 30W purification lamp the luminosity is greater than 300LUX.
Controller	Intelligent button or touch screen
Electric Outlet	Each national standard electric outlet is available
Alt Name	Bench-Top Fume Hood



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com