



TYPE I ULTRAPURE WATER PURIFICATION SYSTEM BKQ-401

TYPE I ULTRAPURE WATER PURIFICATION SYSTEM BKQ-401

PURE / ULTRAPURE WATER SYSTEMS

This series is a high-performance ultrapure water system designed for laboratories. Using tap water as feed, it consistently produces ultrapure water meeting Grade I of GB/T 6682-2008 and CAP requirements.

The system combines two-pass RO technology, a high-rejection RO membrane and a polished stainless-steel sterile tank with built-in UV. Multiple purification stages ensure stable, high-quality water. The 7-inch touchscreen and intelligent controller support online water monitoring and preset-quality, preset-volume dispensing for intuitive and efficient operation.

Its compact split design saves bench space. Quick-change disposable cartridges and dual mechanical/touch dispensing provide convenient and efficient laboratory water access.



1. Compact Design-Easy Operation
7" Touchscreen: Simple control
Split Design: Saves bench space
Dual Dispensing: Mechanical or touch control
2. Two-Pass RO: Sterile Water Tank
High-Rejection RO: Stable performance
Sterile Tank: Stainless steel with built-in UV
3. Reliable Quality: Meets CAP I and GB/T standards
Smart Monitoring: Easy Management
MCU Control: Stores water and quality records
4. Real-Time Display: Water quality and system status
Smart Dispensing: Preset quality and volume

SPECIFICATIONS

Model	BKQ-401
Production Rate	20L/hr
Water Quality	Type I, Type III
Feed Water	Tap Water
Tank Capacity	10L
Working temperature	10-35°C room
Feed Water Quality	Tap water (hardness < 200 ppm)
RO/UP Water Production Capacity	20 L/h
RO/UP Water Dispensing Flow Rate	60 L/h
Pure Water (RO) Production Quality :	
Product Water Conductivity (µS/cm)	≤ 5 µS/cm @ 25°C
Desalination Rate	≥ 99% @ 25°C
Organic Matter Removal Rate	> 99%
Ultrapure Water (UP) Quality :	
Product Water Resistivity (MΩ·cm)	≥ 18.2 MΩ·cm @ 25°C
Total Organic Carbon (TOC)	≤ 2 ppb
Microorganism	≤ 0.1 CFU/mL
Microparticles	≤ 1 particle/mL
Heavy Metal Content	≤ 0.1 ppb
Humidity requirements	30%~85%
Power supply	220V(+10%or-10%),50Hz(+1or-1)Pure(Adapter can be used)
Water tank storage	10 litre
Design Resistivity	18.2MΩ/Cm @25°C
Host size	H 565*W248*D 450mm(Fully drainable UV sterile water tank inside)
Alt Name	Pure / Ultrapure Water Systems

FEATURES

1. Compact Design: Easy Operation
7" Touchscreen: Simple control
Split Design: Saves bench space
Dual Dispensing: Mechanical or touch control
2. Two-Pass RO: Sterile Water Tank
High-Rejection RO: Stable performance
Sterile Tank: Stainless steel with built-in UV
3. Reliable Quality: Meets CAP I and GB/T standards
Smart Monitoring: Easy Management
MCU Control: Stores water and quality records
4. Real-Time Display : Water quality and system status
Smart Dispensing : Preset quality and volume
Stable Performance: Long Service
Constant Flow: Consistent and repeatable output
Quality Piping: Reliable, leak-free operation
5. Easy-to-use operation interface: Equipped with a 7-inch touch human-machine interface, the operation is intuitive and convenient, and users can easily control the operation and parameter settings of the equipment.
6. Efficient water tank design: Use 304 stainless steel mirror sterile water tank, fully sealed method and ultraviolet sterilization technology to ensure the hygiene and purity of water quality.
7. Intelligent control system: The equipment is equipped with an intelligent main control system to achieve fully automatic operation and improve work efficiency and stability.
8. Water quality monitoring guarantee: supports online monitoring of conductivity and resistivity to ensure reliable water quality and avoid using substandard water.
9. Flexible water intake mode: supports qualitative and quantitative water intake to meet different experimental needs and ensure accurate and convenient water intake.

APPLICATIONS

Chromatography & MS: HPLC, UPLC, LC-MS, GC-MS.
Spectroscopy & Elemental Analysis: ICP-MS, ICP-AES, AAS.
Precision Instruments: MALDI-TOF-MS, IC, TOC analyzers.
Chemical Analysis: Electrochemistry, spectrophotometry.
Life Sciences: Media and reagent preparation, especially for low-endotoxin and low-nuclease applications.
General Laboratory Use: Clinical reagent dilution and routine laboratory water needs.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com