



HORIZONTAL LAMINAR AIRFLOW BLHZ-209

HORIZONTAL LAMINAR AIRFLOW BLHZ-209

HORIZONTAL LAMINAR FLOW CABINET



LCD Display, Airflow velocity, UV timer, UV work time, system work time, real time
 Waterproof Socket, 2 waterproof sockets are located in the side panel, for optimum convenience of using small devices inside the cabinet.
 Gas Tap
 UV Lamp, Emission of 253.7 nanometers for most efficient decontamination
 Side Glass Window, The transparent side glass windows maximize light and visibility inside the cabinet, providing a bright and open working environment.

SPECIFICATIONS

Model	BLHZ-209
External Size (W*D*H)	1800*820*2070mm
Internal Size (W*D*H)	1700*500*570mm
Work Surface Height	750mm
Display	LCD Display
Airflow Velocity	Average of 0.3~0.5m/s
Material	Main Body: Cold-rolled steel with anti-bacteria powder coating Work Table: 304 stainless steel Side and Front Window: 5mm toughened glass, anti-UV
Pre-Filter	Polyester fiber, washable
HEPA Filter	99.995% efficiency at 0.3um
Noise	<65dB
Front Window	Motorized
Max Opening	430mm
LED Lamp	16W*2
UV Lamp	30W*1 Emission of 253.7 nanometers
Consumption	500W
Waterproof Socket	Two, total load <500W
Caster	Universal caster with leveling feet
Power Supply	AC220V ±10%, 50/60Hz; 110V±10%, 60Hz
Standard Accessory	LED lamp, UV lamp*2, Base stand, Gas Tap, Waterproof socket*2 Water tap
Optional Accessory	Electric height adjustable base stand
Gross Weight	328kg
Package Size(W*D*H)	1950*1060*1600mm
Alt Name	Horizontal Laminar Flow Cabinet

FEATURES

LCD Display, Airflow velocity, UV timer, UV work time, system work time, real time




Waterproof Socket, 2 waterproof sockets are located in the side panel, for optimum convenience of using small devices inside the cabinet.



Gas Tap



UV Lamp, Emission of 253.7 nanometers for most efficient decontamination



Side Glass Window, The transparent side glass windows maximize light and visibility inside the cabinet, providing a bright and open working environment.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com