



AUTOMATIC SEALING AND CAPPING MACHINE BLIH-101

AUTOMATIC SEALING AND CAPPING MACHINE BLIH-101

CAPMAN AUTOMATIC CAPPING MACHINE

Automatic Sealing and Capping Machine is specifically designed to speed up the process of packaging of your product but also ensure the sterility of the product and that the product does not get contaminated.



Easy to switch the cap opening module, one equipment can achieve 96 holes rack or 100 holes rack cap opening and closing operation.
Quick to switch the cap hat and suitable for all kinds of brand sample tube cap opening and closing.
Built-in integration with an intelligence sensor which can identify different volume sample tubes .
It not only can be independent operation but also can be together with other liquid workstation and robotic arm comprise an automatic system.
Drawer type inlets and outlets for avoiding sample contamination.
Touch screen control panel One-button operation provides the user with direct and convenient operation feeling.

SPECIFICATIONS

Model	BLIH-101
Rack type available	96 or 100 Holes
Screen	Touch Screen control panel
Sensor	Built-in integration with intelligence sensor
Alt Name	Capman Automatic Capping Machine

APPLICATIONS

Pharmaceutical, medical, Food Industry, Beverage Industry.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com