



AUTOMATIC SEALING AND CAPPING MACHINE BLIH-103

AUTOMATIC SEALING AND CAPPING MACHINE BLIH-103

AUTOMATIC SEALER

Automatic Sealer is a full-automatic equipment for sealing PCR plate, deep well plate, micro-plate, etc. It is mainly used for sealed storage of heat-sensitive substances, which can effectively prevent sample loss and cross contamination.



1. Fully automatic and one key operation, no need for additional manual operation.
2. Adapts to varieties of plate types.
3. Adjustable sealing pressure.
4. A variety of communication interfaces are reserved, which can be quickly connected to other machines.
5. Each roll of film can seal more than 350 plates(reference quantity).

SPECIFICATIONS

Model	BLIH-103
Positions	1pc
Method	Cold seal
Speed	<16s/plate
Operating Mode	Manual plate change, automatic sealing
Application	PCR plate, deep well plate and micro-plate conforming to SBS standard (Customizable)
Control Mode	Buttons
Interface	USB
Consumption	75W
Power Supply	100~240V, 50/60Hz
External Size (L*W*H) mm	310*380*410
Package Size (L*W*H) mm	454*382*464
Net Weight	16.5kg
Gross Weight	20.5kg
Alt Name	Automatic Sealer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com