



MULTI-PARAMETER ANALYZER BMET-1109

MULTI-PARAMETER ANALYZER BMET-1109

MULTIMETER

Multi-parameter Analyzer simultaneously measures Ion, Conductivity, Dissolved Oxygen, Temperature that otherwise requires four different meters. With the appropriate probes installed, it can also simultaneously measure several or one parameters.



SPECIFICATIONS

Model	BMET-1109
DC Voltage	
- 6 V / 60 V / 600 V	$\pm (0.8\% + 3)$
AC Voltage	
- 60 mV / 600 mV / 6 V / 60 V / 600 V / 750 V	$\pm (0.8\% + 3)$
Resistance	
- 6 k Ω / 60 k Ω / 600 k Ω	$\pm (1.0\% + 3)$
- 6 M Ω / 10 M Ω	$\pm (1.2\% + 5)$
Capacitance	
60 nF / 600 nF / 6 μ F / 60 μ F / 600 μ F / 6 mF	$\pm (4.0\% + 5)$
Frequency	
30 ~ 1000 Hz	$\pm (0.5\% + 2)$
Display	6000 Counts
Auto Range	Yes
Bargraph Display	Yes
Continuity	Yes
NCV (Non-Contact Voltage)	Yes
Live Test	Yes
Data Hold	Yes
Work Light	Yes
Back Light	Yes
Auto Power Off	Yes
Power Supply	3 x 1.5 V AAA Batteries
Weight	Approx. 243 g
Size	169 x 83 x 53 mm
Safety Standard / Rating	EN61010-1, EN61010-2-033, EN61326, CAT.III 600 V
Alt Name	Multimeter



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com