



PH METER BMET-204

PH METER BMET-204

PORTABLE PH METER

The pH scale typically ranges from 0 to 14. According to the Nernst equation pH values are affected by temperature, hence the temperature compensation is a vital function of any pH Meter to ensure accurate pH measurement.



SPECIFICATIONS

Model	BMET-204
pH	
- Range	-2.00 to 14.00 pH
- Resolution	0.01 pH
- Accuracy	± 0.01 pH
- Calibration Points	Up to 2
- Standard Customization	NIST buffers
mV	
Inputs	
- pH Electrode	BNC(Q9)
Measurement	
- Reading Mode	Continuous
- Reading Prompts	Reading, Stable
- Temp. Compensation	MTC
Display Options	
- Backlight	Yes
- Auto Shutdown	300, 600, 1200, 1800, 3600 sec., off
- IP Rating	IP65
General	
- Power	Rechargeable Lithium batter; AC Adapter, 100-240 V AC input, DC5V output
- Dimensions	80 x 255 x 35 mm
- Weight	400 g (0.88 lb)
Alt Name	Portable pH Meter



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com