



CONDUCTIVITY METER BMET-301

CONDUCTIVITY METER BMET-301

BENCHTOP CONDUCTIVITY METER

A Conductivity meter is an instrument that measures the amount of electric conductivity or current in a solution. The conductance can be measured by applying an alternating electrical current to the two electrodes present in the solution, after which the cations move to the negative electrode and the anions move to the positive electrode. This movement ultimately leads the solution to be conductive.



SPECIFICATIONS

Model	BMET-301
Conductivity	
- Range	0.00 $\mu\text{S}/\text{cm}$ to 100 mS/cm
- Resolution	0.01 $\mu\text{S}/\text{cm}$ minimum; changed with range
- Accuracy	$\pm 1.5\%$ FS
- Reference Temperature	25 $^{\circ}\text{C}$
Measurement	
- Reading Mode	Continuous
- Reading Prompts	Reading, Stable
- Temp. Compensation	MTC
Inputs	
- pH Electrode	BNC(Q9)
- Conductivity Probe	5-pin aviation connector
Display Options	
- Backlight	Yes
- Auto Shutdown	300, 600, 1200, 1800, 3600 sec., off
- IP Rating	IP54
General	
- Power	AC Adapter, 100-240 V AC input, DC9 V output
- Dimensions	200x160x63 mm
- Weight	600 g (1.32 lb)
Alt Name	Benchtop Conductivity Meter



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com