



TURBIDITY METER BMET-701

TURBIDITY METER BMET-701

TURBIDITY METER

Turbidity meters measure the amount of suspended solid particles in water. It works on the backscattering principle, wherein an infrared light is shined into the medium by infrared diode. The electronic device then based on the received signal calculates the relative turbidity of the medium.



SPECIFICATIONS

Model	BMET-701
Measuring Unit	
- Range	(0 ~20.00) NTU, (20.0 ~200.0) NTU
- Resolution	0.01NTU, 0.1NTU, 1NTU
- Accuracy	$\leq \pm 8 \%$
- Fluctuation	$\leq \pm 1.0\% \text{ FS}/30 \text{ min}$
General	
- Power Supply	AA Battery, AC Adapter 200-240 VAC input DC5V output
- Dimension (mm) & Weight (kg)	100x220x80, 0.8
Alt Name	Turbidity Meter



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com