



## MAGNETIC STIRRER BMGS-101

## MAGNETIC STIRRER BMGS-101

### MINI MAGNETIC STIRRER

Mini Magnetic Stirrer is a mini-sized and portable magnetic stirrer with high standard reliability and durability. It has wide range of application.



Distinct designed. Small and exquisite. Simple operation.  
Long using life with no carbon brush motor.  
Speed regulation stepless, large range of speed.  
Strong magnetism to avoid stirrer jumping.  
Convenient and portable. Meet different kinds of experimental demands.

## SPECIFICATIONS

Model	BMGS-101
Platform Size	Φ88mm
Platform Material	Plastic
Speed Range	0~2500 rpm
Max. Capacity	250ml
Stirring bar Dimension	Φ8x46mm
Motor Input/Output	1.2W/0.8W
Voltage	DC12V 1A
Power	≤12W
Dimension	W.110xD.110xH.70mm
Net Weight	1.2kg
Alt Name	Mini Magnetic Stirrer



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)