



## MAGNETIC STIRRER BMGS-106

# MAGNETIC STIRRER BMGS-106

MAGNETIC STIRRER



Digital speed control within the range of 0-1500 rpm

Maximum stirring quantity of H<sub>2</sub>O at 3L

LED display shows stirring speed

## SPECIFICATIONS

Model	BMGS-106
Work plate Dimension	Ø135mm (5 inch)
Work plate material	Nylon+GF
Motor type	Brushless DC motor
Stirring positions	1
Max. stirring quantity [H <sub>2</sub> O]	3L
Max. magnetic bar [length]	50mm
Speed range	200-1500rpm
Speed display	LED
Protection class	IP42
Motor output Power	3W
Power	15W
Voltage, Frequency	100-240V, 50/60Hz
Dimension [WxDxH]	150x260x80mm
Weight	1.8kg
Permissible ambient temperature and humidity	10-40°C, 80%RH
Alt Name	LED Digital Magnetic Stirrer



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)