



## MAGNETIC STIRRER BMGS-108

# MAGNETIC STIRRER BMGS-108

MAGNETIC STIRRER



Speed range of 0-1500 rpm  
Maximum stirring quantity of H<sub>2</sub>O up to 10L  
Ceramic work plate provides excellent corrosion resistance

## SPECIFICATIONS

Model	BMGS-108
Work plate Dimension [W x D]	184x184mm (7 inch)
Work plate material	Glass ceramic
Motor type	Shaded pole motor
Stirring positions	1
Max. stirring quantity [H <sub>2</sub> O]	10L
Max. magnetic bar [length]	80mm
Speed range	0-1500rpm
Speed display	Scale
Protection class	IP21
Motor output Power	1.5W
Power	30W
Voltage, Frequency	100-120/200-240V, 50/60Hz
Dimension [WxDxH]	215x360x112mm
Weight	3.8kg
Permissible ambient temperature and humidity	5-40°C, 80%RH
Alt Name	7" Square Magnetic Stirrer



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)