



## MAGNETIC STIRRER BMGS-109

# MAGNETIC STIRRER BMGS-109

ULTRA-FLAT COMPACT MAGNETIC STIRRER



Speed range 15 to 1500 rpm  
Automatic reverse rotation for excellent mixing  
PET top work plate provides excellent chemical resistance  
Motorless technology provides ultra silent and vibration-free operation  
Perfect magnetic stirrer for use during titration (bottle-top burette) in analytical chemistry

## SPECIFICATIONS

Model	BMGS-109
Work plate Dimension	Ø90mm
Work plate material	PET
Stirring positions	1
Max. stirring quantity[H <sub>2</sub> O]	0.8L
Max.magnetic bar[length]	40mm
Speed range	15-1500rpm
Protection class	IP21
Output Power	3W
Power	5W
Voltage, Frequency	100-240V, 50/60Hz
Dimension[WxDxH]	175x120x15mm
Weight	0.3kg
Permissible ambient temperature and humidity	5-40°C, 80%RH
Alt Name	Ultra-flat Compact Magnetic Stirrer



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)