



BIOLOGICAL MICROSCOPE BMIC-601

BIOLOGICAL MICROSCOPE BMIC-601

MULTI-FUNCTION DIGITAL BIOLOGICAL MICROSCOPE

The Biological Microscope is a powerful scientific instrument for studying the microscopic world of biology. It offers high-resolution imaging, versatile observation modes, and precision controls, making it an indispensable tool for researchers, educators, and enthusiasts.



SPECIFICATIONS

Model	BMIC-601
Optical System	/
Viewing Head	Siedentopf binocular viewing head, inclined at 30°,
Eyepiece	WF10X/18
Objective	Achromatic objective 4X,10X, 40X,100X
Nosepiece	Quadruple nosepiece
Focusing	Coaxial coarse&fine adjustment, fine division 0.002mm, coarse stroke 37.7mm per rotation, fine stroke 0.2mm per rotation, moving range 20mm
Stage	Double layers mechanical stage: 132*142mm, moving range: 75*40mm
Condenser	Abbe condenser NA1.25 with aperture diaphragm
Illumination	S-LED Illumination, brightness adjustable
Power Supply	110~220V, 50/60Hz
Accessories	Photo attachment, video attachment
Package Size	323*264*440mm
Gross Weight	6.5kg
Alt Name	Multi-function Digital Biological Microscope



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com