



METALLURGICAL MICROSCOPE BMIC-803

METALLURGICAL MICROSCOPE BMIC-803

METALLURGICAL MICROSCOPE

Metallographic microscope is a high-tech product which combines optical microscope technology, photoelectric conversion technology and computer image processing technology perfectly. It can easily observe metallographic images on a computer, analyze and grade metallographic maps, edit, save and print the pictures.



SPECIFICATIONS

Model	BMIC-803
Viewing Head	Sliding Trinocular Head Inclined at 45°
Nosepiece	Quadruple Nosepiece
Eyepiece	WF10X18mm
Objective	Achromatic Objectives 4X, 10X, 40X, 100X
Storage	Double Layer Mechanical Stage 140*135mm moving range 70*50mm
Focusing	Coaxial coarse&fine focusing adjustment with rack and pinion mechanism fine focusing scale value 0.002mm
Illumination	Build-in Illumination, halogen lamp 6V/20W
Filter	Blue, Green, Yellow
Power Supply	AC110~220V, 50/60Hz
Package Size(mm)	330*280*440 / 370*260*470
Gross Weight(kg)	6.5; 5
Alt Name	Metallurgical Microscope



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com