



## MICROPLATE READER BMRW-102

# MICROPLATE READER BMRW-102

## MICROPLATE READER

Microplate Reader is a laboratory instrument that is used to measure chemical, biological or physical reactions, properties and analytes within the well of a microplate. Flexibility and reliability makes this product an unique choice.



Bichromatic, 8-channel Optical system.  
Large storage capacity up to 200 test protocols and 100,000 results.  
LED light source with lifespan over 10,000 hours.  
3 shaking modes with variable time.  
Windows operation system, self-check, malfunction alarm.  
7 inch touch TFT LCD screen, test results simultaneously displayed in screen.

## SPECIFICATIONS

Model	BMRW-102
Repeatability	<1%
Stability	<0.005 Abs/10 min
Linearity	R $\geq$ 0.995
Display	7 inch TFT LCD, 800x480 Pixels
Interface	RS~232, 4 USB ports, SD Card
Printer	Built-in printer, external printer
Power Supply	AC 100~240V, 50/60 Hz
Dimension	400Lx310Wx169H mm
Weight	6 kg
Readout Range	0.000~3.500 Abs
Plate Type	96-well plate
Light Source	Monochromatic LED
Reading Speed	5s with single wavelength for whole plate
Input	Mouse, membrane keyboard, touch screen
Humidity	15%~85%
Temperature	10°C~35°C
Measure System	8~channel optical system
Band Width and Accuracy	8 nm $\pm$ 2 nm
Shaking	3 shaking modes with variable time
Wavelength	405, 450, 492, 630nm, two more optional
Alt Name	Microplate Reader

## APPLICATIONS

Drug Discovery, Research, Laboratory, Medical, Assay Development, Life Science.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)