



## MULTISTATION STIRRER BMSS-101

## MULTISTATION STIRRER BMSS-101

### MULTI-POSITION MAGNETIC STIRRER

Multi-position Magnetic Stirrer is designed for the high-efficiency experiments and has 4/8/12 stirring bar configurations to choose. It can use laboratory space efficiently. DC brushless motor makes the model, low speed but stable, high speed but powerful with low noise. Stainless steel is coated by a silicone resin film disk. Non-slip and anti-corrosion.



Digital display stirring speed at real-time, stirring bars operate synchronously and provide uniform agitation reaction conditions.

Efficient multi-position magnetic stirrers configure 4/8/12 stirring positions to use lab space efficiently.

DC brushless motor. Low speed but stable, high speed but powerful. Low noise and free maintenance.

Stainless steel is coated by silicone resin film pad. Non-slip and anti-corrosion.

## SPECIFICATIONS

Model	BMSS-101
Stirring Point Quantity	4
Size of Stirring Bars	Φ8x46mm
Max. Stir Capacity	400mlx4
Digital Display	Speed Display
Speed Range	200~1200rpm
Temp. Control Range (surface)	/
Platform Material	stainless steel with silicon pad
Motor	DC brushless motor
Voltage	AC100~230V, 50/60Hz
Power	20W
Fuse	250V,1A, Φ5x20
Dimension	W.115xD.480xH.50mm
Net Weight	3kgs
Alt Name	Multi-Position Magnetic Stirrer



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)