



## MULTISTATION STIRRER BMSS-104

# MULTISTATION STIRRER BMSS-104

## MAGNETIC STIRRER



Efficient, simple and safe stirring  
High performance and space-saving design to maximize the use of benchtop space and saving time by stirring multiple solutions  
Maintenance-free brushless DC motor  
Speed range 0-1100 rpm  
Stainless steel work plate covered with silicone cushion, provides excellent performance of skid resistance

## SPECIFICATIONS

Model	BMSS-104
Work plate Dimension [W x D]	180x450mm
Work plate material	Stainless steel with silicone coating
Motor type	Brushless DC motor
Stirring positions	10
Max. stirring quantity[H <sub>2</sub> O]	0.4Lx10
Max.magnetic bar[length]	40mm
Speed range	0-1100rpm
Speed display	Scale
Protection class	IP20
Motor output Power	4W
Power	20W
Voltage, Frequency	100-120/200-240V, 50/60Hz
Dimension[WxDxH]	622x205x65mm
Weight	3.2kg
Permissible ambient temperature and humidity	5-40°C, 80%RH
Alt Name	Magnetic Stirrer



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)